



angiodroid

ANGIODROID
Carbon Dioxide
Interventional Peripheral
Angiography



Contents

- [1- Angiodroid: The CO² Injector](#)
- [2 - Angiodroid: business opportunities](#)
- [3 - Interventional Peripheral Angiography](#)
- [4 - Carbon Dioxide: Gas Properties](#)
- [5 - CO2 Indications & Applications](#)
- [6 - How it works](#)
- [7 - The Market: Overview](#)
- [8 - Clinical validation and scientific resources](#)
- [9 - Carbon Dioxide Toxicity Analysis](#)
- [10 - Certification](#)
- [11 - Brochure](#)
- [12 - CO2 Imaging](#)
- [13 - Carbon Dioxide Procedures - Video](#)
- [14 - Carbon Dioxide History](#)
- [15 - Angiodroid: conclusion](#)



Angiodroid: The CO² Injector



angiodroid

Is a ***CO₂ Contrast Media Injector***
dedicated to
***Interventional and Diagnostic
Angiographic Peripheral Procedures***



Angiodroid: Main Fields of Application



angiodroid

has 3 main Fields of Application

- **Interventional Radiology**
- **Vascular Surgery**
- **Interventional Cardiology**



Angiodroid: conclusion



angiodroid

Carbon Dioxide (CO²) is the only Safe Contrast Agent for patient with

- **Hyper sensibility to Iodinated Contrast Material**
- **Renal Failure**
- **Diabetic Diseases**



Angiodroid: conclusion



angiodroid

Is the only Device Worldwide on Market

- **100% DIGITAL**
- **100% AUTOMATIC**
- **100% SAFE**
- **100% USER FIENDLY**



Angiodroid: Business Opportunities

Single-use Kit
connection line +
4 Frenchs Catheter



Angiodroid The Injector



After sale services
- Periodical maintenance
- CO2 cylinder substitution



Interventional Peripheral Angiography

- Rising life Expectancy
- Increment of vascular diseases
- Development of sophisticated radiological interventional techniques

**All those Aspects have determined an increment
of Angiographic Procedures**



CIN – Contrast Induced Nephropathy

Contrast-Induced Nephropathy (CIN)

Each angiographic procedure requires the injection of iodinated contrast media that, in an increasing number of procedures, has aspects of clear contraindications or quantitative limitation: the case of diseases associated with impaired renal or diabetic vascular dysfunctions.

Contrast-Induced Nephropathy (CIN) is a form of Acute Renal Failure caused by exposure to contrast media during cardiology and radiology cathetering procedures. The lack of effective treatment to prevent CIN remains problematic for patients with renal-insufficiency.



CIN – Contrast Induced Nephropathy

What are the potential consequences?

CIN is associated with increased:

- Major adverse in-hospital cardiac events
- In-hospital mortality rates
- Long-term mortality
- Risk of acceleration toward end-stage renal disease (dialysis)
- Longer and more frequent hospital stays

According to “Marenzi, et al”, hospitalized patients who received contrast media and who acquired CIN had significantly higher mortality rate (31% vs. 0.6%) than patients who did not acquire CIN.

CIN risk incidence: 30%  **higher health care costs**




CO₂ as Solution to CIN Risk

An alternative to the use of iodinated contrast media is represented by Carbon Dioxide (CO₂), a gas injected into the vessels results in a clear contrast effect and is easily removed from the circulation via the lungs.


In the '70s, Hawkins was the pioneer in the US of CO₂ angiography for high risk patients with intolerance to iodine or renal failure.

With the advent of image subtraction angiography in 1980, the use of low-density CO₂ became possible, and then, through the implementation of digital technologies for image subtraction, angiography with CO₂ became an applicable method in the different fields of angiography.




CO2 and Iodine angiography venograms comparison

CO2



a.

Iodine



b.

Figure 1: Anteroposterior venograms of the elbow region in a 27-year-old man with renal insufficiency and a failed hemodialysis AVF in the contralateral wrist and elbow. (a) CO₂ venogram shows the cephalic (white arrows), basilic (arrowheads), and deep (black arrows) veins. (b) Conventional venogram depicts only the cephalic (arrows) and basilic (arrowheads) veins.



Carbon Dioxide: Gas Properties

Gas Properties:

- Highly compressible
- 20 times more soluble than O²
- Non Viscous and Buoyant
- Invisible, Colourless and Odourless
- Radiopacity
- Rapidly dissolving in Blood
- Lacks both Allergic Potential & Renal Toxicity

Next:

- CO² doesn't mix with Blood
- It has extreme diffusibility
- None Total Maximum dosage



CO2 angiography: critical issues

Basically there are 3 critical issues in the use of CO2 as a contrast medium:

1. to determine the amount of injected gas and to remove the air from all over the system, in order to avoid the risk of vascular emboli
2. to inject the gas at controlled pressure in order to avoid vessel rupture in the presence of obstructions (aneurysms)
3. to properly modulate the gas injection; in order to ensure emptying of the catheter from the internal fluid, at an early stage, and to avoid the "jet", while controlling the gas injection in a second stage

The technological solutions available up to today (see *The Market: Overview*), do not solve these critical issues by relying on manual control of the operator or neglecting the variability of blood pressure during the injection phase. In fact the amount of gas injected through a syringe of known volume depends the pressure inside the vessel which is not constant over time.



CO₂ Angiography Indication

Indications

- renal failure diabetes
- intolerance to iodinated contrast
- creatinine greater than 1.8 mg / dl

Contraindications

- above diaphragm arterial studies
- pulmonary insufficiency, pulmonary AV malformation
- interauricular or interventricular communication

Possible side effects

- nausea
- pain
- dizziness
- tachycardia



CO₂ Clinical Application

Diagnosis applications:

- pelvic arteries
- lower extremity
- renal, visceral, venous studies

Interventional applications:

- angioplasty
- stent placement
- shunt for hemodialysis
- vena cava filter placement
- ablation of renal artery
- transcatheter embolization
- endovascular treatment of abdominal aortic aneurysm



How it works

Injection Technique and Imaging

1. Use of Closed Circuit
2. Pressure and Volume Control System
3. Selective Injection
4. Use of Nitro-glycerine (0,1 mg) for Blood low range Patients
5. Wasteless Diet to eliminate interferences
6. Administer Buscopan (20 mg) to reduce bowel movement
7. 3-4 frames/sec Imaging Acquisition and 60 ms
8. Use Fluoroscopy mode
9. Wait 1-2 min between different injections
10. Use “stacking” software during Post processing
11. Trendelenburg modality while delivering (10° - 15°)




How it works

Combined Approach

CO2


Indication



- ✓ Borderline Renal Failure
- ✓ Chronic Nephropathy
- ✓ Allergic Potential
- ✓ Riduzione carico rene



Iodinated Contrast Agent

Indication




- ✗ Thoracic Angiography Procedures
- ✗ Coronary Angiography
- ✗ Neurovascular Studies
- ✗ Non cooperative patients
- ✗ Rear vessels

Peripheral Angiography
PTCA
Endoprostesisys



THE VERSATILE CONTRAST & FLUID MANAGEMENT SYSTEM



ANGIOFLUSH III fluid collection system

- manual injection kit
- operator-dependent: no procedural repeatability



CO2 ANGIOSET

- manual injection kit
- difficulty of achieving the desired injection volume





CADDI

- automatic volume adjustment
- no control of injection pressure
- withdrawn from the market because of its complex usability
- none technical innovations: it's a simply mechanical handling of a syringe

* Recalled Device



Angiodroid: The CO² Injector



angiodroid

100% DIGITAL

100% AUTOMATIC

100% SAFETY

100% USER FIENDLY



ANGIODROID: advantages and uniqueness

Major Technical figures:

- Definition of the amount of gas to be injected into the vessel
- Definition of Pressures and Injection volumes (automation, security and repeatability of the procedure)
- Total removal of the air contained in the catheter through the initial injection of a small quantity of gas
- Constant injection pressure of gas during the procedure, removing the risk of vessel rupture in presence of obstructive aneurysms
- Positive pressure of internal pneumatic circuit reducing to zero the risk of introducing air



Injection Systems Comparison

Manual Injection

- Anti-bacterial Filter
- Syringe
- Valves
- Check Valves
- Special Connectors

HANDICAP

- ✘ High X-ray Operator Dose Exposure
- ✘ Operator dependent
- ✘ Serious contamination risk with ambient air
- ✘ Difficult to accurately determine the injected dose
- ✘ Increased risk of pain for the patient

Digital Automatic Injection

Digital Automatic injector that accurately determines Volume and Pressure


UNIQUENESS

- ✓ Reduced X-ray exposure
- ✓ Locked and Safe System
- ✓ Automatic cleaning of Internal circuit
- ✓ User friendly
- ✓ None risks of Air Ambient Contamination



Manual Injection

The CO₂ manual Injection doesn't allow an optimal Gas release control, it exposes to the risk of Contamination and requires a long training period




Pressure/time graph of two manual CO₂ injections. Phase A corresponds to line washing; phase B corresponds to gas injection –Corazza et al 2013²



Digital Automatic Injection

Angiodroid The CO₂ Injector has exceeded all those problems (Angiodroid®, Angiodroid SRL, Bologna, ITALIA)




P1 and P2 pressure recordings during automated line washout and CO₂ gas injection. (Corazza et. al 2013)



Digital Automatic Injection

The amount of injected gas should be proportional to the volume of the vessel to have in view, while the gas pressure, which regulates the gas flow, should be proportional to blood flow in the vessel.

- Vessel volume
- Gas volume
- Blood flow
- Gas pressure



Clinical validation and scientific resources

NCBI Resources How To

PubMed.gov
US National Library of Medicine
National Institutes of Health

PubMed carbon dioxide angiography Search

RSS Save search Limits Advanced

Display Settings: Summary, 20 per page, Sorted by Recently Added

Send to:

Results: 1 to 20 of 721

<< First < Prev Page 1 of 37 Next > Last >>

[Safety of carbon dioxide digital subtraction angiography.](#)

1. Moos JM, Ham SW, Han SM, Lew WK, Hua HT, Hood DB, Rowe VL, Weaver FA.
Arch Surg. 2011 Dec;146(12):1428-32.

PMID: 22288088 [PubMed - indexed for MEDLINE]

[Related citations](#)

[Rising PaCO₂ in the ICU: using a physiologic approach to avoid cognitive biases.](#)

2. Schwartzstein RM, Parker MJ.
Chest. 2011 Dec;140(6):1638-42. No abstract available.

PMID: 22147823 [PubMed - indexed for MEDLINE]

[Related citations](#)

[Mechanical aspects of CO₂ angiography.](#)

3. Corazza I, Rossi PL, Feliciani G, Pisani L, Zannoli S, Zannoli R.
Phys Med. 2011 Dec 2. [Epub ahead of print]

PMID: 22138139 [PubMed - as supplied by publisher]

[Related citations](#)

[Discontinuation of the plastic bag delivery system for carbon dioxide angiography will increase radiocontrast nephropathy and life-threatening complications.](#)

4. Cho KJ, Hawkins IF Jr.
AJR Am J Roentgenol. 2011 Nov;197(5):W940-1. No abstract available.

PMID: 22021546 [PubMed - indexed for MEDLINE]

[Related citations](#)



Clinical Validation

Physica Medica

In Press, Corrected Proof - Note to users

doi:10.1018/j.ejmp.2011.11.003 | How to Cite or Link Using DOI

[Permissions & Reprints](#)



ALMA MATER STUDIORUM
UNIVERSITÀ DI BOLOGNA

Original Paper

Mechanical aspects of CO₂ angiography

Ivan Corazza^a, Pier Luca Rossi^a, Giacomo Feliciani^a, Luca Pisani^b, Sebastiano Zannoli^b, Romano Zannoli^a

^a Cardiovascular Dept., University of Bologna, Via Massarenti, 9, 40138 Bologna, Italy

^b Spark Srl, Bologna, Italy

Received 5 July 2011; revised 13 October 2011; Accepted 12 November 2011. Available online 3 December 2011.

Abstract

The aim of this paper is to clarify some physical–mechanical aspects involved in the carbon dioxide angiography procedure (CO₂ angiography), with a particular attention to a possible damage of the vascular wall.

CO₂ angiography is widely used on patients with iodine intolerance. The injection of a gaseous element, in most cases manually performed, requires a long training period. Automatic systems allow better control of the injection and the study of the mechanical behaviour of the gas.

CO₂ injections have been studied by using manual and automatic systems. Pressures, flows and jet shapes have been monitored by using a cardiovascular mock. Photographic images of liquid and gaseous jet have been recorded in different conditions, and the vascular pressure rises during injection have been monitored.

The shape of the liquid jet during the catheter washing phase is straight in the catheter direction and there is no jet during gas injection. Gas bubbles are suddenly formed at the catheter's hole and move upwards: buoyancy is the only governing phenomenon and no bubbles fragmentation is detected. The pressure rise in the vessel depends on the injection pressure and volume and in some cases of manual injection it may double the basal vascular pressure values.

CO₂ angiography is a powerful and safe procedure which diffusion will certainly increase, although some aspects related to gas injection and chamber filling are not yet well known. The use of an automatic system permits better results, shorter training period and limitation of vascular wall damage risk.

Keywords: Angiography; CO₂ angiography; Interventional radiology



Carbon Dioxide Toxicity Analysis

Dogs accomplished Studies

Renal Study

7cc/kg of CO₂ were injected into the renal artery in lateral decubitus

- **Radioisotope infusion**, to determine renal function, was done:
 - 24 hours pre-injection
 - Immediately after injection
 - 24 hours after injection
- Immediate Post CO₂ injection - Transient change in blood flow
- 24 hours Post CO₂ injection - Resumed normal flow and function
- SEM showed no changes to endothelial lining, renal parenchyma

Dr. Hawkins, Univ. of Florida

Conclusion: CO₂ has no renal toxicity

Carbon Dioxide Toxicity Analysis

Rabbits accomplished Studies

Hepatic Study

CO₂ was injected into the hepatic artery

- Blood samples to determine ALT (Alanine Transferase)
 - Pre injection (Baseline – Normal)
 - Immediately post –injection (Transient elevation)
 - 7 days post-injection (Return to baseline normal)
- Liver pathology study 7 day post-injection:
 - No significant abnormalities in histopathology

Dr. Mladinich, Univ. of Florida

Conclusion: CO₂ has no long term hepatic toxicity.

Carbon Dioxide Toxicity Analysis

Pigs accomplished Studies

Venous Study

15 CO₂ / 5 IC, (0.2cc/kg - 6.4cc/kg) IVC

- CO₂ injected in three different positions :
 - Supine
 - Right decubitus
 - Left decubitus
- Parameters obtained pre-injection and 1, 3, 5 , 10 minutes after injection:
 - Pulmonary arterial pressure *Transient increase*
 - Systemic arterial pressure *No change*
 - Arterial blood gases(1,3,5,10min. Post inj.) *No change*
- 2 Expired (1IC 3.2cc/kg, 1CO₂ 6.4cc/kg rd)
- ***Conclusion: Diagnostic doses have no significant effect on cardiopulmonary function.***

Dr. Cho - Univ. of Michigan

Carbon Dioxide Toxicity Analysis

Dogs accomplished Studies

Coronary/Neuro Study

4-9 injections for a total of 120-300cc per dog

- 9 - injections ascending aorta
- 3 - injections common carotid
- 2 - injections both arteries

Shifrin, et al

No Change in EEG, ECG, Arterial Blood Gases.

Carbon Dioxide Toxicity Analysis

Rabbits accomplished Studies

Coronary / Neuro Study

27 received 12cc/kg CO₂

11 received 12cc/kg Renografin 76

11 received 12cc/kg Saline

All injections into the left ventricle

Dr. Bettmann Boston Univ

- ✓ ***No damage to Myocardium***
- ✓ ***No damage to Blood/Brain barrier***
- ✓ ***No differences among study groups***

Carbon Dioxide Toxicity Analysis

Rats accomplished Studies

Intra-carotid injection of CO₂ produces:

- Multifocal ischemic infarctions
 - Disrupted blood-brain barrier
 - Lesions of the endothelial cell membrane
- *Conclusion: CO₂ should not be used for angiographies of cerebral arteries*

Carbon Dioxide Toxicity Analysis

Human accomplished Studies

Retrospective Study

208 procedures / 189 patients (138 males / 70 female, 9 to 86 yrs.)

- CO₂ injected in all 208 procedures
- CO₂ and iodinated contrast 175 procedures
- CO₂ alone 33 procedures.

Results:

- 32 AE with CO₂, 6 with iodinated contrast
- Most common CO₂ AE: GI related (nausea)
- Most common with iodinated AE: Allergic reaction

University of Florida

Carbon Dioxide Toxicity Analysis

Human accomplished Studies

Safety Study

21 Patients, ages of 50 to 90 years

10 patients CO₂ / 11 patients Optiray® 320.

Safety: AE's, ABG levels of CO₂, O₂, pH, bicarb levels

Efficacy: paired images read by investigator

Results:

- ✓ 1 AE with iodinated contrast, none with CO₂
- ✓ No significant changes (pH, CO₂, O₂, or bicarbonate)
- ✓ All CO₂ films diagnostic, (IC Films - 45% diagnostic-55% excellent)

Boston University and Dartmouth

Conclusion: CO₂ is at least as Safe as Iodinated Contrast Medium.

Carbon Dioxide Toxicity Analysis

Retrospective Study

Human accomplished

- ❑ **1997 – 2007 - 8 years**

- ❑ **Individual Operator**

- ❑ **654 procedures**
 - **245 angiographies**
 - **64 CO2 alone**
 - **181 Combined approach (CO2 e Iodine)**
 - **409 PTCA**
 - **185 CO2 alone**
 - **224 Combined approach (CO2 e Iodine)**

Carbon Dioxide Toxicity Analysis

Retrospective Study

Human accomplished

Results

- ❑ **119 pain and annoyance attacks within 7290 injections**
 - ❑ 75 on the first 180 injections per a total of 12 procedures
 - ❑ 44 on the followind 7110 injection, within 642 procedures

- ❑ **482 100% positive results (No IODINE)**
 - ❑ Alone CO2 use
 - ❑ Relevant Imaging

- ❑ **157 Partial positive results (IODINE < 50 ml)**
 - ❑ Good Imaging Quality of some Arterial Region
 - ❑ Combined approach required (with IODINE)

- ❑ **15 failure results (IODINE > 50 ml)**
 - ❑ Non cooperative Patients
 - ❑ Exagerate Interference with Bowel movement
 - ❑ Low Imaging Quality

Carbon Dioxide Toxicity Analysis

Carbon Dioxide and Iodine Combines Study

- ❑ 82 patients
- ❑ Random Prospective Study
- ❑ Renal Angiograph and PTR
- ❑ Combined Approach vs IODIO alone

The quantity of delivered Iodinated Contrast Agent is strictly related to Creatinine increase within 2 days post procedure. (**p=0.011**)

Higher is the use of Iodine, Higher is CIN Risk

Procedures	Creatinine before Procedures	Creatinine After 3 days	Creatinine Significance
IODINE 100% failure	1,8 ± 0,4	2,9 ± 0,9	P < 0.05
IODINE partial failure	1,9 ± 0,6	2.0 ± 0,8	ns
NO IODINE 100% success	1,7 ± 0,8	1,8 ± 0,4	ns

Carbon Dioxide Toxicity Analysis

Carbon Dioxide Angiography Effective Analysis

❑ 50 cases – Retrospective Study

PTA-PTRA

❑ 88% with CO₂

❑ 12% combined approach with a small Iodine quantity *Kessel DO et al: Cardiovasc Interv Rad 25(6):476-83 2002*

❑ Renal Transplantation Analysis

❑ 17 patients – Retrospective Study

❑ None false information

❑ 100% correct diagnosis

Chao, Major, Weaver et al. :J Vascular Surgery 45(3);451-60 2007

❑ 100 cases . Retrospective Study

EVAR

❑ CO₂ vs IODIO

❑ None Creatinine change

❑ Most part of successfully procedures

Lorsch et al.: CO₂ angio in transplant Rontgenpraxis 55(1) 26-32 2003

Carbon Dioxide Toxicity Analysis

Animal accomplished studies

- ✓ No Renal Toxicity
- ✓ No Epatic Toxicity
- ✓ No significant effect on Cardiopulmonary function
- ✗ Possible neurotoxicity on rats. None evidence on dogs and rubbitts

Human accomplished studies (more than 300 pubblications)

- ✓ Retrospective study: 208 procedures/ 189 pat → 32 (16%) AE (Bowel movement)
- ✓ Safety evaluation: 61 pat → none changements
ph, CO₂, O₂, HCO₃



CE mark: medical device

Notified body: IMQ

Applicant: SIAS Spa

CERTIFICATO CE
Certificato n. 1497/MDD

Dichiarazione di approvazione del sistema qualità
(Sistema completo di garanzia qualità)

Visto l'esito delle verifiche condotte in conformità all'Allegato II della direttiva 93/42/CEE e s.m.i., si dichiara che la ditta:

S.I.A.S. S.r.l.
40057 CADRIANO DI GRANAROLO EMILIA (BO) - VIA MINGHETTI 9/11 (ITA) - Italy

mantiene negli stabilimenti di:

40050 FUNO DI ARGELATO (BO) - VIA ZANOTTI 4 (ITA) - Italy

un sistema qualità che assicura la conformità dei seguenti prodotti:

Iniettore per mezzo di contrasto
Mod. ANGIODROID
Marca SIAS

ai requisiti essenziali della direttiva suddetta ad essi applicabili (in tutte le fasi dalla progettazione al controllo finale).

Riferimento pratiche IMQ: 10AL00102.

Questa Dichiarazione di approvazione è rilasciata dall'IMQ S.p.A. quale organismo notificato per la direttiva 93/42/CEE e s.m.i.
Il numero identificativo dell'IMQ S.p.A. quale organismo notificato è: 0051.

Emesso il: 2012-01-15


IMQ

Questa Dichiarazione di approvazione è soggetta alle condizioni previste da IMQ nel "Regolamento per la certificazione CE dei dispositivi medici in base alla direttiva 93/42/CEE".
Essa non è comunque valida dopo il 2017-01-15 (articolo 11, comma 11 della direttiva).

 IMQ

IMQ S.p.A. Società a ruolo unico - 20136 Milano - Via Quintano 43 - tel. (02) 7511131 - fax (02) 5091500 - info@imq.it - www.imq.it
REA MI 1209284 - Registro Imprese MI 1209284/0010 - C.F./P.I. 12092840100 - capitale sociale € 100.000 euro



CE mark: single-use connection line

Notified body: TUV

Applicant: Bioengineering Laboratories Spa

AY12/06 ZERTIFIKAT ◆ CERTIFICATE ◆ 認証証書 ◆ CERTIFICADO ◆ CERTIFICAT


Product Service

CERTIFICATO CE
Garanzia di qualità della produzione
Direttiva 93/42/CEE concernente i dispositivi medici (DDM), allegato V
(dispositivi in classe IIa, IIb o III)
N° G2 12 02 44964 028

Fabbricante: Bioengineering Laboratories S.p.A.
Via Isonzo 9/B
20036 Meda (MI)
ITALIA

Categoria(e) di prodotto: Cateteri per gastromanometria, rettali, uretrali, urologici e loro accessori; ago per blocco paracervicale uterino, circuiti per l'afèresi del sangue; Clamp vascolari; linee di estensione per angiografia

Con il presente certificato, l'organismo notificato di TÜV SÜD Product Service GmbH certifica che il fabbricante sopra menzionato ha implementato un sistema di qualità per la fabbricazione ed il controllo finale dei dispositivi / categorie di dispositivi in questione secondo quanto stabilito nella direttiva DDM, allegato V. Questo sistema di qualità risponde ai requisiti della presente direttiva ed è soggetto a sorveglianze regolari. Per l'ammissione sul mercato di dispositivi delle classi IIb e III, è richiesto un certificato addizionale, di cui all'allegato III. Osservare le note riportate sul retro.

N° del rapporto: ITA 2071251S

Valido da: 2012-02-21
Valido fino al: 2016-03-02


Data: 2012-02-22
Hans-Heiner Junker


306807

TÜV SÜD PRODUCT SERVICE GMBH è Organismo Notificato con numero identificativo 0123.
Pagina 1 di 2
Traduzione per scopi informativi. La sola versione inglese è legalmente impegnativa.

TÜV SÜD Product Service GmbH · Zertifizierungsstelle · Ridlerstraße 55 · 80339 München · Germany 



**Practical applications of the advanced Angiodroid
analgesic, nontoxic CO₂ contrast media system**

CO₂ is used in various types of peripheral vascular interventions, such as:

- Angioplasty
- Stent placement
- Transcatheter embolization
- Treatment of abdominal aortic endovascular aneurysm
- Renal ablation

CO₂ is also used as a contrast media in cases of:

- Cholangiography
- Nephrostomy
- Gastrostomy
- Visualization of abscesses and cysts

In the vascular field CO₂ is the ideal contrast media in the following cases:

- Pelvic arteries
- Lower leg arteries
- Renal arteries
- Visceral Arteries
- Visualization of arteriovenous fistulas for hemodialysis
- Venous Studies
- Retrograde Portography

References:

1 Radiology Volume 211 Number 1 - October 2008
2 Journal of the American College of Cardiology
Volume 46 Issue 5 - November 2005 Pages 752-758
Luh Shih-wei MD, Gary D. Webb MD, Petera M. Kelly MD,
George Sengen MD, FRC, Resonance J. Lundy MD, et al



The Gada Group Limited, Grosvenor Gardens House, 35-47, Grosvenor Gardens,
SW1W 0BS London, Phone: +44 0 20 7802 0011 • Fax: +44 0 20 7828 1211
info@gadagroup.com gadagroup.com



angiodroid



Angiodroid Innovative CO₂ System

**Resolving the critical limitations of peripheral angiography
with clearly visible benefits for patients**

Angiodroid carbon dioxide peripheral interventional angiography: Eliminating critical issues for improved patient safety

● Peripheral interventional angiography overview.


The rise in life expectancy, the increase in vascular disease and the development of increasingly sophisticated interventional radiological techniques is leading to a progressive increase in investigative procedures in cardiac angiography (coronary angiography and PTCA) in the brain and peripheral areas.

Each angiographic procedure requires the injection of iodinated contrast medium which features clear contraindications or quantitative limitations: especially in the case of conditions associated with impaired renal function or diabetic vascular disease. Patients who present these chronic conditions (more than 20% of cases treated) are particularly susceptible to the toxic effects of contrast agents and to developing Contrast-Induced Nephropathy (CIN) - a form of Acute Renal Failure.

It is, therefore, necessary to find alternative methods that address the risk of intolerance for the patient and reduce the economic impact to the medical community.

● CO₂: the easier, safer and established alternative system.

An alternative, 'CIN-FREE', media is represented by carbon dioxide (CO₂), a gas injected into the vessels that results in a clear contrast effect and is easily removed from the circulation via the lungs.



Anteroposterior venograms of the elbow region in a 27-year-old man with renal insufficiency and a failed hemodialysis AVF in the contralateral wrist and elbow.
(a) CO₂ venogram shows the cephalic (white arrows), basilic (black arrows), and deep (black arrows) veins.
(b) Conventional venogram depicts only the cephalic (white arrow) and basilic (black arrows) veins.*

● Angiodroid is the new, patented system for CO₂ injection that resolves the critical limitations associated to most devices currently on the market.

- Angiodroid is 100% safe and eliminates the risk of CIN.
- Angiodroid allows the interventionalist to define and precisely control the amount of gas to be injected into the vessel.
- Angiodroid guarantees complete removal of the air contained in the injection catheter, providing total patient safety.
- Angiodroid ensures a constant injection pressure of the gas during the procedure with extreme ease and precision, allowing the repeatability of the procedure to obtain quality images while eliminating the risk of vessel rupture in the presence of obstructive aneurysms.
- Angiodroid uses a very thin catheter (CO₂ is about 400 times less viscous than iodine) allowing more selective catheterization.
- Over 100 clinical reviews and publications demonstrate the effectiveness of CO₂ angiography.



Example of image quality as seen in modern based digital technologies.

● Compatible with most modern digital imaging systems

The Angiodroid injection technique has been optimized for use with modern, currently available, angiographic systems that implement algorithms for digital image acquisition.



CO2 Imaging



CO2 Imaging




CO2 Imaging




CO2 and Iodine angiography venograms comparison

CO2



a.

Iodine



b.

Figure 1: Anteroposterior venograms of the elbow region in a 27-year-old man with renal insufficiency and a failed hemodialysis AVF in the contralateral wrist and elbow. (a) CO₂ venogram shows the cephalic (white arrows), basilic (arrowheads), and deep (black arrows) veins. (b) Conventional venogram depicts only the cephalic (arrows) and basilic (arrowheads) veins.



Carbon Dioxide History

CO² has been in use in Radiology as an imaging agent since 1914 to evaluate the retroperitoneum and to image kidneys and tumors.

CO² was in the routinely in the early 1950s for detection of pericardial effusion.

In the 70s Dr. Hawkins and Dr. Cho pioneered the intraarterial use of CO² . With the advent of Digital Subtraction Angiography (DSA) in 1980, reliable imaging of "low density" contrast material became available. With the addition of high-resolution DSA, Stacking SW, Tilting Tables, and reliable delivery systems, CO² is used not only in patients with contrast material allergies and renal failure, but also in patients undergoing routine angiographic studies and many interventional procedures



Angiodroid: conclusion



Carbon Dioxide (CO²) is the only **Safe Contrast Agent** for patient with

- ***Hypersensitivity to Iodinated Contrast Material***
- ***Renal Failure***
- ***Diabetic Diseases***

Carbon Dioxide can be injected as Contrast Agent in any Luminal Structure (Peripheral Regions mandatory) except for thoracic aorta and its branches



Angiodroid: conclusion



angiodroid

Is the only Device Worldwide on Market

- 100% DIGITAL
- 100% AUTOMATIC
- 100% SAFETY
- 100% USER FIENDLY





angiodroid

ANGIODROID
Carbon Dioxide
Interventional Peripheral
Angiography

Thanks for the attention

